beeem
elevate | everything | energetically

EDUCATIONAL SERIES

Tattoos & Scars

Understanding the intended energetic approach

Why tattoos and scars are discussed in energetic wellness

Tattooing repeatedly punctures the skin and deposits pigment beneath its surface. Scars form as skin repairs after an injury or procedure. Tattoo inks vary in composition, and some can cause allergic reactions or other complications. These physical issues are distinct from the energetic concepts described below. [1]

Some holistic wellness models describe tattoos and scars as possible sources of energetic interference or “static” within the body’s communication pathways. A kink in a garden hose is a useful analogy for this proposed model; it is not a demonstrated biological mechanism. A healed tattoo or scar does not, by itself, establish that an energetic block exists.

Step 1 | Neutralize with beeem ON™

beeem describes ON™ as structured green light carrying frequency information inspired by natural materials. Its intended energetic role is to help neutralize a perceived stressor associated with a tattoo or scar and support energetic coherence, without aiming to change the tattoo or scar physically.

Step 2 | Elevate with beeem IN™

After ON™, beeem’s protocol uses IN™ structured red light, described as carrying a distinct frequency library inspired by botanicals, minerals, crystals, and other natural materials. This second step is intended to introduce supportive energetic information. These are descriptions of frequency information, not the delivery of botanical ingredients, minerals, or antioxidants into the skin.

What has been reported?

In practitioner reports supplied by beeem, tattoos described as energetic stressors during Quantum Reflex Analysis (QRA) no longer tested that way after 60 seconds of ON™ exposure. These preliminary muscle-response observations do not establish toxin removal, restored biological pathways, clinical benefit, or permanence. They also do not establish the added benefit of IN™ or equivalent results for scars.

beeem
elevate | everything | energetically

EDUCATIONAL SERIES

Tattoos & Scars

ON™ first. IN™ second. 60 seconds each.

The beeem™ application protocol

The sequence below reflects beeem's supplied protocol for existing, fully healed tattoos and scars. The 60-second duration is a brand protocol, not a clinically validated treatment dose.

1 Neutralize | Green beeem ON™

Direct the green light over the entire tattoo or scar for 60 seconds. For a larger area, divide it into sections and apply ON™ to each section separately for 60 seconds.

2 Elevate | Red beeem IN™

Immediately follow with the red light over the same area for an additional 60 seconds. For a larger area, apply IN™ to each section separately for 60 seconds.

3 The area has been beeemed™

The two-step sequence is complete: first ON™ for intended energetic neutralization, then IN™ for intended energetic elevation. beeem presents this as a one-time application for an existing tattoo or scar; a lasting effect has not been established by the reports provided. New tattoos or new areas of trauma are considered separately after healing.

What this process does - and does not - mean

The intended focus is the energetic experience associated with the area. beeem™ is not a tattoo-removal or scar-revision procedure and is not intended to alter ink, color, scar tissue, or appearance. "Neutralize" does not mean sterilize the skin, remove metals or chemicals, or eliminate the medical risks of tattoo ink.

Care and responsible use

Follow the device's safety instructions and eye-protection requirements; never direct the light toward eyes. Do not use this guide as wound-care advice for fresh tattoos, open skin, or recent surgery. Continue professional aftercare, and seek medical advice for persistent pain, redness, swelling, rash, or other concerns.

These statements describe the intended energetic use of beeem™ technology and are not intended to diagnose, treat, cure, or prevent any disease. This educational guide does not replace appropriate medical care. Results and experiences may vary.

[1] Physical tattoo background: FDA, Tattoos & Permanent Makeup: Fact Sheet. This source does not evaluate beeem™. Product descriptions, timing, and QRA observations: beeem-supplied information; no controlled study data provided for this guide.